

mdi 0.1µm AseptiCap WS-y are low protein binding hydrophilic PVDF membrane capsule filters, validated to retain mycoplasma, a critical requirement for sterilization of mammalian cell culture media.

These capsules offer serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane for enhanced throughput.

0.1µm AseptiCap WS-y capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Very low hold up volume in filters

Applications

Sterile Filtration of

- Cell culture media
- Growth regulators
- Small Volume Parenterals

Ordering Information




Microbially Validated as per ASTM F 838-05
Complies with USFDA 21 CFR 210.3(b)(6)
Meets and Exceeds USFDA 21 CFR 177.1520

Specifications

Materials of Construction

Membrane	Hydrophilic PVDF
Support Layer	Polyester
Plastic Components	Polypropylene

Typical Water Flow Rates



Microbial Retention

LRV >7 for *Acholeplasma laidlawii* (ATCC 23206) per cm²

Maximum Operating Temperature

80 °C @ ≤30 psi (2 Kg/cm²)

Maximum Differential Pressure

60 psi (4 Kg/cm²) @ 30 °C

Bubble Point with 50% IPA

≥31psi

Sterilization

By Irradiation: Gamma irradiatable up to 50 kGy

By Autoclave: Autoclavable at 125°C for 30 minutes, 1 cycle after gamma irradiation. Can not be in-line steam sterilized

Toxicity

Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics

Cytotoxicity

Passes Biological Reactivity Tests, *In Vitro*, USP <87> for cytotoxicity

Bacterial Endotoxin

Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>

Fiber Release

Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release

Particle Release

The filtrate complies with USP <788> test for particulate matter in injections

TOC and Conductivity

Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after flushing with a specified volume of WFI

Extractables with WFI

Passes NVR test as per USP <661>

Oxidizable Substances

Passes test as per USP <1231>

Type	Size		Pore Size	Inlet /Outlet		Radiation Sterilizable		Inline/T-Line		Sterility		Pack Size	
	Code	Length and EFA	Code	Code	Code	Yes	Code	Inline	Code	Code	Code	1	Code
AseptiCap WS (0.2µm upstream)	LWS1	5" (3000 cm ²)	53	0.1µm	36	1 1/2" Triclover	E	Yes	R	X	Non-Sterile	1	01
		10" (6000 cm ²)	54			Single Step 1/2" Hose Barb	Q	No**	X	T-Line***	T	Gamma Sterile	3
AseptiCap WS (0.45µm upstream)	LWSX	20" (12000 cm ²)	55			3/8" Hose Barb	I						
		30" (18000 cm ²)	56			1" Hose Barb*	Z						

*1" Hose Barb is not available in 5" capsule filters

**Gamma sterile capsule filters cannot be gamma irradiated again

***T-line Capsule Filter are available with 1 1/2" Sanitary Flange I/O Connections only

Note: Size 5" is available in Inline Capsule filters only

EFA: Effective Filtration Area

Example

LWS1	54	36	QQ	R	X	1	01
------	----	----	----	---	---	---	----

For Non-Sterile: LWS15436QQRX101

For Gamma Sterile: LWS15436QQXX301

DST LWSX36R1630B

mdi 0.1 µm AseptiCap WS-γ are low protein binding hydrophilic PVDF membrane capsule filters, validated to retain mycoplasma, a critical requirement for sterilization of mammalian cell culture media.

These capsules offer serial filtration incorporating a larger pore size upstream membrane to protect the downstream membrane for enhanced throughput.

0.1 µm AseptiCap WS-γ capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Very low hold up volume in filters

Applications

Sterile Filtration of

- Cell culture media
- Growth regulators
- Small Volume Parenterals

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size				
	Code	Length and EFA	Code		Code		Code		Code			Code		Code			
AseptiCap WS (with 0.45µm Upstream)	DWSX	1" (100 cm²)	31	0.1µm	36	¼" SHB	A	Yes	R		Non-Sterile	1	1	01			
		1" (250 cm²)	51			½" Hose Barb	D	No****	X		Gamma Sterile	3					
AseptiCap WS (with 0.2µm Upstream)	DWS1	2" (500 cm²)	52			1½" Sanitary Flange	E	* Single step ½" hose barb and 3/8" hose barb end connections are not available in 1" Capsule filter **Male luer slip end connection is available as outlet only in 1" capsule filters ***3/16" hose barb end connection is available in: - 1" and 2" capsule filters as inlet and outlet - 5" as outlet only ****Gamma sterile capsule filters cannot be gamma irradiated again									
		5" (1000 cm²)	53			¾" Sanitary Flange	S										
8" (2000 cm²)		57	Quick Connector			J											
Single Step ½" Hose Barb*		Q															
						Female Luer Lock	U										
						Male Luer Slip**	W										
						3/16" Hose Barb***	N										
						3/8" Hose Barb*	I										
EFA: Effective Filtration Area																	
Example																	
DWSX		53		36		QQ		R		X	1		01				

For Non-Sterile: DWSX5336QQRX101

For Gamma Sterile: DWSX5336QQXX301

DST DWSX36R1620C



Specifications

Materials of Construction

Membrane	Hydrophilic PVDF
Plastic Components	Polypropylene

Microbial Retention


LRV >7 for *Acholeplasma laidlawii* (ATCC 23206) per cm²

Bubble Point with 50% IPA/Water
≥31psi

Maximum Operating Temperature
80 °C @ ≤30 psi (2 Kg/cm²)

Maximum Differential Pressure
60 psi (4 Kg/cm²) @ 30 °C

Typical Water Flow Rates



Sterilization

By Irradiation: Gamma irradiatable up to 50 kGy

By Autoclave: Autoclavable at 125°C for 30 minutes, 1 cycle after gamma irradiation. Can not be in-line steam sterilized

Toxicity

Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics

Cytotoxicity

Passes Biological Reactivity Tests, In Vitro, USP <87> for cytotoxicity

Bacterial Endotoxin

Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>

Fiber Release

Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release

Particle Release:

The filtrate complies with USP <788> test for particulate matter in injections

TOC and Conductivity

Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after a 3 liter flush

Extractables with WFI

Passes NVR test as per USP <661>

Oxidizable Substances

Passes test as per USP <1231>

Complies with USFDA 21 CFR 210.3(b)(6)

Meets and Exceeds USFDA 21 CFR 177.1520

mdi AseptiCap WS are low protein binding hydrophilic PVDF membrane inline capsule filters, designed for sterile filtration of very small fluid volumes in formulation and process development labs.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.



Applications

Sterile Filtration of

- Cell culture media
- Growth regulators
- Small Volume Parenterals

Complies with USFDA 21 CFR 210.3(b)(6)

Meets and Exceeds USFDA 21 CFR 177.1520

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Very low hold up volume

Specifications

Construction		
Pore Size	0.1µm	
Membrane	Hydrophilic PVDF	
Plastic Components	Polypropylene	
Size		
Size	25 mm	50 mm
Effective Filtration Area (Nominal)	5 cm²	20 cm²
Integrity Testing/Retention		
Bubble Point	≥ 31psi (2.18Kg/cm²) with 50% IPA/Water Solution	
Microbial Retention	LRV>7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm²	
Operational		
Max. Operating Temperature	55 °C	60 °C
Max. Differential Pressure	75 psi (5 Kg/cm² @25°C)	42 psi (3 Kg/cm²) @ 30 °C
Sterilization	By Gas	Sterilization by Ethylene Oxide
	By Autoclave	Autoclavable at 125°C for 30 minutes, 2 cycles. Cannot be in-line steam sterilized
Assurance		
Bacterial Endotoxin	Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>	
Toxicity	Passes Biological reactivity Test, <i>In Vivo</i> , as per USP <88> for Class VI plastics	

Assurance

Cytotoxicity	Passes Biological Reactivity Tests, <i>In Vitro</i> , USP <87> for cytotoxicity
Fiber Release	Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release
Particle Release	The filtrate complies with USP <788> test for particulate matter in injections
TOC and Conductivity	Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after flushing with a specified volume of WFI
Extractables with WFI	Passes test as per USP <661>
Oxidizable Substances	Passes test as per USP <1231>

Ordering Information

25 mm Inline Capsule Filters

Type		Size		Pore Size		Inlet		Outlet		X	X	Sterility		Pack Size	
	Code		Code		Code		Code		Code				Code		Code
<i>AseptiCap WS</i> (0.45µm Upstream)	IWSX	25mm	06	0.1 µm	36	1/8" Hose Barb	H	1/8" Hose Barb	H			Non Sterile	1	100	04
<i>AseptiCap WS</i> (0.2µm Upstream)	IWS1					1/4" Hose Barb	B	1/4" Hose Barb	B			EO Sterile	2		
						Female Luer Lock	M	Male Luer Slip	N						
								Male Luer Lock	L						

Example:

IWSX	06	36	M	N	X	X	1	04
------	----	----	---	---	---	---	---	----

50 mm Inline Capsule Filters

Type		Size		Pore Size		Inlet		Outlet		X	X	Sterility		Pack Size	
	Code		Code		Code		Code		Code				Code		Code
<i>AseptiCap WS</i> (0.45 µm Upstream)	IWSX	50 mm	10	0.1 µm	36	1/4" SHB	B	1/4" SHB	B			Non Sterile	1	10	02
<i>AseptiCap WS</i> (0.2 µm Upstream)	IWS1					3/4" Sanitary Flange	S	3/4" Sanitary Flange*	S			EO Sterile	2	100	04
<i>Vented AseptiCap WS</i> (0.45 µm Upstream)	VWSX														
<i>Vented AseptiCap WS</i> (0.2 µm Upstream)	VWS1														

* In vented *AseptiCap WS* 3/4" Sanitary Flange is available as outlet only

Example:

IWSX	10	36	S	S	X	X	1	04
------	----	----	---	---	---	---	---	----

MDI AseptiCap WS are low protein binding hydrophilic PVDF membrane capsule filters offering serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane for enhanced throughput.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Low hold up volume

Applications

- Cell culture media
- Growth regulators
- Small volume parenterals

Complies with USFDA 21 CFR 210.3(b)(6)

Meets and Exceeds USFDA 21 CFR 177.1520



Specifications

Pore Size

0.1 µm

Materials of Construction

Membrane	Hydrophilic PVDF
Support Layer	Polyester
Plastic Components	Polypropylene

Microbial Retention

LRV >7 for *A. laidlawii* (ATCC 23206) per cm²

Maximum Operating Temperature

80 °C @ ≤ 30 psi (2 Kg/cm²)

Maximum Differential Pressure

60 psi (4 Kg/cm²) @ 30 °C

Bubble Point

≥31psi (2.18 Kg/cm²) with 50% IPA/ water solution

Sterilization

By Autoclave: Autoclavable at 125°C for 30 minutes, 2 cycles. Can not be in-line steam sterilized

By Gas: Sterilization by Ethylene Oxide

Toxicity

Passes Bioreactivity test, *In Vivo*, as per USP <88> for Class VI plastics

Cytotoxicity

Passes Biological Reactivity Tests, *In Vitro*, USP <87> for cytotoxicity

Bacterial Endotoxin

Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>

Fiber Release

Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release

Particle Release:

The filtrate complies with USP <788> test for particulate matter in injections

TOC and Conductivity

Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after flushing with a specified volume of WFI

Extractables with WFI

Passes test as per USP <661>

Oxidizable Substances

Passes test as per USP <1231>

Ordering Information

Type	Size	Pore Size	Inlet /Outlet	X	Inline/T-Line	Sterility	Pack Size
Code	Length and EFA	Code	Code		Code	Code	Code
AseptiCap WS (0.2 µm upstream)	5" (3000 cm ²)	53	0.1 µm	36	1½" Sanitary Flange	E	
	10" (6000 cm ²)	54			Single Step ½" Hose Barb	Q	
AseptiCap WS (0.45 µm upstream)	20" (12000 cm ²)	55			3/8" Hose Barb	I	
	30" (18000 cm ²)	56			1" Hose Barb	Z	

EFA: Effective Filtration Area

*T-line Capsule Filter are available with 1½" Sanitary Flange I/O Connections only
Size 5" is available in Inline Capsule filters only

Example

LWS1	56	36	QQ	X	X	1	01
------	----	----	----	---	---	---	----

DST LWSX36X1600B

Advanced Microdevices Pvt. Ltd., 20-21, Industrial Area, Ambala Cantt - 133006, INDIA

Tel: +91-171-2699290, 2699471 Email: info@mdimembrane.com Website: www.mdimembrane.com

mdi AseptiCap WS are low protein binding hydrophilic PVDF membrane capsule filters offering serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane for enhanced throughput.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Very low hold up volume in filters

Applications

Sterile Filtration of

- Cell culture media
- Growth Regularors
- Small Volume Parenterals

Complies with USFDA 21 CFR 210.3(b)(6)

Meets and Exceeds USFDA 21 CFR 177.1520



Specifications

Pore Size

0.1 μ m

Materials of Construction

Membrane	Hydrophilic PVDF
Support Layer	Polyester
Plastic Components	Polypropylene

Microbial Retention

LRV>7 for *A.laidlawii* (ATCC 23206) per cm^2

Maximum Operating Temperature

80 °C @ \leq 30 psi (2 Kg/ cm^2)

Maximum Differential Pressure

60 psi (4 Kg/ cm^2) @ 30 °C

Bubble Point

\geq 31 psi (2.18Kg/ cm^2) with 50% IPA/water Solution

Sterilization

By Autoclave: Autoclavable at 125°C for 30 minutes, 2 Cycles. Can not be in-line steam sterilized

By Gas: Sterilization by Ethylene Oxide

Toxicity

Passes Bioreactivity test, *In Vivo*, as per USP <88> for Class VI plastics

Cytotoxicity

Passes Biological Reactivity Tests, *In Vitro*, USP <87> for cytotoxicity

Bacterial Endotoxin

Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>

Fiber Release

Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release

Particle Release:

The filtrate complies with USP <788> test for particulate matter in injections

TOC and Conductivity

Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after a 3 liter flush

Extractables with WFI

Passes NVR test as per USP <661>

Oxidizable Substances

Passes test as per USP <1231>

Ordering Information

Type		Size		Pore Size		Inlet /Outlet		X	X	Sterility		Pack Size	
	Code	Length and EFA	Code		Code		Code				Code		Code
AseptiCap WS (0.2 μ m Upstream)	DWS1	1" (250 cm^2)	51	0.1 μ m	36	1/4" SHB	A			Non-Sterile	1	1	01
		2" (500 cm^2)	52			1/2" Hose Barb	D			EO Sterile	2		
AseptiCap WS (0.45 μ m Upstream)	DWSX	5" (1000 cm^2)	53			1 1/2" Sanitary Flange	E						
		8" (2000 cm^2)	57			3/4" Sanitary Flange	S						
						Quick Connector	J						
						Single Step 1/2" Hose Barb*	Q						
						Female Luer Lock	U						
						Male Luer Slip**	W						
						3/16" Hose Barb***	N						
						3/8" Hose Barb*	I						
Example		DWSX	53		36	EE		X	X		1		01

*Single step 1/2" hose barb and 3/8" Hose Barb end connections are not available in 1" capsule filter

**Male Luer slip end connection is available as outlet only in 1" capsule filters

***3/16" hose barb end connection is available in:
- 1" and 2" capsule filters as inlet and outlet
- 5" as outlet only

mdi AseptiCap WS- γ are low protein binding hydrophilic PVDF gamma sterilizable membrane inline capsule filters, designed for sterile filtration of very small fluid volumes in formulation and process development labs.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.



Applications

Sterile Filtration of

- Cell culture media
- Growth regulators
- Small Volume Parenterals

Complies with USFDA 21 CFR 210.3(b)(6)

Meets and Exceeds USFDA 21 CFR 177.1520

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Very low hold up volume

Specifications

Construction		
Pore Size	0.1 μm	
Membrane	Hydrophilic PVDF	
Plastic Components	Polypropylene	
Size		
Size	25 mm	50 mm
Microbial Retention	LRV>7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm²	
Integrity Testing/Retention		
Bubble Point (with 50% IPA/Water)	≥ 31 psi (2.18Kg/cm²)	
Effective Filtration Area (Nominal)	5 cm²	20 cm²
Operational		
Max. Operating Temperature	55 °C	60 °C
Max. Differential Pressure	75 psi (5 Kg/cm² @25°C)	42 psi (3 Kg/cm²) @ 30 °C
Sterilization	By Irradiation	Gamma Irradiatable up to 50 kGy
	By Autoclave	Autoclavable at 125°C for 30 minutes, 1 cycle after gamma irradiation. Cannot be in-line steam sterilized
Assurance		
Bacterial Endotoxin	Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test	
Toxicity	Passes Biological reactivity Test, <i>In Vivo</i> , as per USP <88> for Class VI plastics	
Cytotoxicity	Passes Biological Reactivity Tests, <i>In Vitro</i> , USP <87> for cytotoxicity	

Assurance

Fiber Release	Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release
Particle Release	The filtrate complies with USP <788> test for particulate matter in injections
TOC and Conductivity	Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after flushing with a specified volume of WFI
Extractables with WFI	Passes test as per USP
Oxidizable Substances	Within limits as specified in USP

Ordering Information

25 mm Inline Capsule Filters

Type		Size		Pore Size		Inlet		Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code			Code		Code
AseptiCap WS (0.45µm Upstream)	IWSX	25mm	06	0.1 µm	36	½" Hose Barb	H	½" Hose Barb	H	Yes	R		Non Sterile	1	100	04
AseptiCap WS (0.2µm Upstream)	IWS1					¼" Hose Barb	B	¼" Hose Barb	B	No	X		Gamma Sterile	3		
						Female Luer Lock	M	Male Luer Slip	N							
								Male Luer Lock	L							

Example:

IWSX	06	36	M	N	X	X	1	04
------	----	----	---	---	---	---	---	----

Note: Gamma Sterile filters can not be sterilized again

Example for Non Sterile: IWSX0636MNRX104

Example for Gamma Sterile: IWSX0636MNXX304

50 mm Inline Capsule Filters

Type		Size		Pore Size		Inlet		Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code			Code		Code
AseptiCap WS (0.45 µm Upstream)	IWSX	50 mm	10	0.1 µm	36	¼" SHB	B	¼" SHB	B	Yes	R		Non Sterile	1	10	02
AseptiCap WS (0.2 µm Upstream)	IWS1					¾" Sanitary Flange	S	¾" Sanitary Flange*	S	No	X		Gamma Sterile	3	100	04
Vented AseptiCap WS (0.45 µm Upstream)	VWSX															
Vented AseptiCap WS (0.2 µm Upstream)	VWS1															

* In vented AseptiCap WS ¾" Sanitary Flange is available as outlet only

Example:

IWSX	10	36	SS	SS	R	X	1	04
------	----	----	----	----	---	---	---	----

Note: Gamma Sterile filters can not be sterilized again

Example for Non Sterile: IWSX1036SSRX104

Example for Gamma Sterile: IWSX1036SSXX304

mdi AseptiCap WS-γ are low protein binding hydrophilic PVDF gamma sterilizable membrane capsule filters offering serial filtration incorporating a larger pore size upstream membrane to protect the downstream membrane for enhanced throughput.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Low hold up volume

Applications

Sterile Filtration of

- Antibodies
- Protein Solutions
- Buffers
- Vaccine concentrates
- Small Volume Parenterals

Microbially Validated as per ASTM F 838-05

Complies with USFDA 21 CFR 210.3(b)(6)

Meets and Exceeds USFDA 21 CFR 177.1520



Specifications


Pore Size

0.2 μm and 0.45 μm

Materials of Construction

Membrane	Hydrophilic PVDF
Plastic Components	Polypropylene

Typical Water Flow Rates, 0.2μm Filters



Microbial Retention

0.2μm: LRV >7 for *B. diminuta* (ATCC 19146) per cm²

0.45μm: LRV >7 for *Serratia marcescens* (ATCC 14756) per cm²

Maximum Operating Temperature

80 °C @ ≤30 psi (2 Kg/cm²)

Maximum Differential Pressure

60 psi (4 Kg/cm²) @ 30 °C

Bubble Point

0.2 μm: > 50psi (3.51Kg/cm²) with Water

0.45 μm: > 30 psi (2.11 Kg/cm²) with water

Sterilization

By Irradiation: Gamma irradiatable up to 50kGy

By Autoclave: Autoclavable at 125°C for 30 minutes, 1 cycle after gamma irradiation. Cannot be in-line steam sterilized

Toxicity

Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics

Cytotoxicity

Passes Biological Reactivity Tests, *In Vitro*, USP <87> for cytotoxicity

Bacterial Endotoxin

Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>

Fiber Release

Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release

Particle Release

The filtrate complies with USP <788> test for particulate matter in injections

TOC and Conductivity

Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after flushing with a specified volume of WFI

Extractables with WFI

Passes test as per USP <661>

Oxidizable Substances

Within limits as specified in USP <1231>

Ordering Information

Type	Size		Pore Size		Inlet /Outlet		Radiation Sterilizable		Inline/T-Line		Sterility		Pack Size		
	Code	Length and EFA	Code		Code		Code		Code		Code		Code		Code
AseptiCap WS (0.45 µm upstream)	LWSX	5" (3000 cm ²)	53	0.2µm	01	1½" Triclover	E	Yes	R	Inline	X	Non-Sterile	1	1	01
		10" (6000 cm ²)	54	0.45µm	02	Single Step ½" Hose Barb	Q	No**	X	T-Line***	T	Gamma Sterile	3		
		20" (12000 cm ²)	55			3/8" Hose Barb	I	**0.45 µm capsule filters are available with 0.8 µm upstream only **Gamma sterile capsule filters cannot be gamma irradiated again ***T-line Capsule Filters are available with 1½" Sanitary Flange I/O Connections only Size 5" is available in Inline Capsule filters only							
AseptiCap WS (0.8 µm upstream)*	LWSS	30" (18000 cm ²)	56			1" Hose Barb	Z								
EFA: Effective Filtration Area															

**0.45 μm capsule filters are available with 0.8 μm upstream only
 **Gamma sterile capsule filters cannot be gamma irradiated again
 ***T-line Capsule Filters are available with 1½" Sanitary Flange I/O Connections only
 Size 5" is available in Inline Capsule filters only

Example

LWSX	56	01	QQ	R	X	1	01
------	----	----	----	---	---	---	----

For Non-Sterile: LWSX5601QQRX101

For Gamma Sterile: LWSX5601QQXX301

DST LWSXRX1600B

mdi AseptiCap WS are low protein binding hydrophilic PVDF membrane inline capsule filters, designed for sterile filtration of very small fluid volumes in formulation and process development labs.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.



Microbially Validated as per ASTM F 838-05
Complies with USFDA 21 CFR 210.3(b)(6)
Meets and Exceeds USFDA 21 CFR 177.1520

Applications

Sterile Filtration of

- Antibodies
- Protein Solutions
- Buffers
- Heat labile additives
- Vaccine concentrates
- Small Volume Parenterals

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Very low hold up volume

Specifications

Construction		
Pore Size	0.2µm	
Membrane	Hydrophilic PVDF	
Plastic Components	Polypropylene	
Size		
Size	25 mm	50 mm
Effective Filtration Area (Nominal)	5 cm²	20 cm²
Integrity Testing/Retention		
Bubble Point	≥ 50psi (3.52 Kg/cm²) with Water	
Microbial Retention	LRV>7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm²	
Operational		
Max. Operating Temperature	55 °C	60 °C
Max. Differential Pressure	75 psi (5 Kg/cm² @25°C)	42 psi (3 Kg/cm²) @ 30 °C
Sterilization	By Gas	Sterilization by Ethylene Oxide
	By Autoclave	Autoclavable at 125°C for 30 minutes, 2 cycles. Cannot be in-line steam sterilized
Assurance		
Bacterial Endotoxin	Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>	
Toxicity	Passes Biological reactivity Test, <i>In Vivo</i> , as per USP <88> for Class VI plastics	

Assurance

Cytotoxicity	Passes Biological Reactivity Tests, <i>In Vitro</i> , USP <87> for cytotoxicity
Fiber Release	Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release
Particle Release	The filtrate complies with USP <788> test for particulate matter in injections
TOC and Conductivity	Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after flushing with a specified volume of WFI
Extractables with WFI	Passes test as per USP <661>
Oxidizable Substances	Passes test as per USP <1231>

Ordering Information

25 mm Inline Capsule Filters

Type		Size		Pore Size		Inlet		Outlet		X	X	Sterility		Pack Size	
	Code		Code		Code		Code		Code				Code		Code
AsepticCap WS (0.45µm Upstream)	IWSX	25mm	06	0.2 µm	01	1/8" Hose Barb	H	1/8" Hose Barb	H			Non Sterile	1	100	04
						1/4" Hose Barb	B	1/4" Hose Barb	B			EO Sterile	2		
						Female Luer Lock	M	Male Luer Slip	N						
								Male Luer Lock	L						

Example:

IWSX	06	01	M	N	X	X	1	04
------	----	----	---	---	---	---	---	----

50 mm Inline Capsule Filters

Type		Size		Pore Size		Inlet		Outlet		X	X	Sterility		Pack Size	
	Code		Code		Code		Code		Code				Code		Code
AsepticCap WS (0.45 µm Upstream)	IWSX	50 mm	10	0.2 µm	01	1/4" SHB	B	1/4" SHB	B			Non Sterile	1	10	02
						3/4" Sanitary Flange	S	3/4" Sanitary Flange*	S			EO Sterile	2	100	04
Vented AsepticCap WS (0.45 µm Upstream)	VWSX														

* In vented AsepticCap WS 3/4" Sanitary Flange is available as outlet only

Example:

IWSX	10	01	S	S	X	X	1	04
------	----	----	---	---	---	---	---	----

mdi AseptiCap WS-γ are low protein binding hydrophilic PVDF gamma sterilizable membrane capsule filters offering serial filtration incorporating a larger pore size upstream membrane to protect the downstream membrane for enhanced throughput.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Very low hold up volume in filters

Applications

Sterile Filtration of

- Antibodies
- Protein Solutions
- Buffers
- Vaccine concentrates
- Small Volume Parenterals

Microbially Validated as per ASTM F 838-05

Complies with USFDA 21 CFR 210.3(b)(6)

Meets and Exceeds USFDA 21 CFR 177.1520

Ordering Information

Type	Size	Pore Size	Inlet /Outlet	Radiation Sterilizable	X	Sterility	Pack Size
Code	Length and EFA	Code	Code	Code		Code	Code
AseptiCap WS (with 0.45µm Upstream)	1" (100 cm ²)	31	0.2µm 01	¼" SHB	A	Yes	R
	1" (250 cm ²)	51	0.45µm 02	½" Hose Barb	D	No****	X
AseptiCap WS (with 0.8µm Upstream)*	2" (500 cm ²)	52	1½" Triclover	E	*0.8 µm upstream is available with 0.45µm capsule filters only **Single step ½" hose barb and 3/8" Hose Barb end connections are not available in 1" capsule filter ***Male luer slip end connection is available as outlet only in 1" capsule filters ****3/16" hose barb end connection is available in: - 1" and 2" capsule filters as inlet and outlet - 5" as outlet only *****Gamma sterile capsule filters cannot be gamma irradiated again		
	5" (1000 cm ²)	53	¾" Sanitary Flange	S			
	8" (2000 cm ²)	57	Quick Connector	J			
			Single Step ½" Hose Barb**	Q			
			Female Luer Lock	U			
			Male Luer Slip***	W			
			3/16" Hose Barb****	N			
			3/8" Hose Barb**	I			
DWSX	53	01	QQ	R	X	1	01

EFA: Effective Filtration Area

Example

For Non-Sterile: DWSX5301QQRX101

For Gamma Sterile: DWSX5301QQXX301




Specifications

Pore Size
0.2 µm

Materials of Construction

Membrane	Hydrophilic PVDF
Plastic Components	Polypropylene

Typical Water Flow Rates, 0.2µm



Microbial Retention

0.2 µm: LRV >7 for *B. diminuta* (ATCC 19146) per cm²

0.45µm: LRV >7 for *Serratia marcescens* (ATCC 14756) per cm²

Maximum Operating Temperature

80 °C @ ≤30 psi (2 Kg/cm²)

Maximum Differential Pressure

60 psi (4 Kg/cm²) @ 30 °C

Bubble Point (with water)

0.2 µm: ≥ 50psi (3.51 Kg/cm²)

0.45 µm: ≥ 30 psi (2.11 Kg/cm²)

Sterilization

By Irradiation: Gamma irradiatable up to 50 kGy

By Autoclave: Autoclavable at 125°C for 30 minutes, 1 Cycle after gamma irradiation. Cannot be in-line steam sterilized

Toxicity

Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics

Cytotoxicity

Passes Biological Reactivity Tests, In Vitro, USP <87> for cytotoxicity

Bacterial Endotoxin

Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>

Fiber Release

Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release

Particle Release

The filtrate complies with USP <788> test for particulate matter in injections

TOC and Conductivity

Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after flushing with a specified volume of WFI

Extractables with WFI

Passes NVR test as per USP <661>

Oxidizable Substances

Passes test as per USP <1231>

mdi AseptiCap WS are low protein binding hydrophilic PVDF membrane capsule filters offering serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane for enhanced throughput.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

Key features

- Absolute retention
- 100% integrity tested
- Low protein binding
- Low extractables
- Very low hold up volume in filters

Applications

Sterile Filtration of

- Antibodies
- Protein Solutions
- Buffers
- Vaccine concentrates
- Small Volume Parenterals

Microbially Validated as per ASTM F 838-05

Complies with USFDA 21 CFR 210.3(b)(6)

Meets and Exceeds USFDA 21 CFR 177.1520



Specifications


Pore Size

0.2 µm and 0.45 µm

Materials of Construction

Membrane	Hydrophilic PVDF
Support Layer	Polyester
Plastic Components	Polypropylene

Typical Water Flow Rates, 0.2 µm



Microbial Retention

0.2 µm: LRV >7 for *B. diminuta* (ATCC 19146) per cm²

0.45 µm: LRV >7 for *S. marcescens* (ATCC 14756) per cm²

Maximum Operating Temperature

80°C @ ≤30 psi (2 Kg/cm²)

Maximum Differential Pressure

60 psi (4 Kg/cm²) @ 30°C

Bubble Point (with water)

0.2 µm: ≥ 50psi (3.51Kg/cm²)

0.45 µm: ≥ 30 psi (2.11 Kg/cm²)

Sterilization

By Autoclave: Autoclavable at 125°C for 30 minutes, 2 Cycles. Can not be in-line steam sterilized

By Gas: Sterilization by Ethylene Oxide

Toxicity

Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics

Cytotoxicity

Passes Biological Reactivity Tests, In Vitro, USP <87> for cytotoxicity

Bacterial Endotoxin

Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>

Fiber Release

Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release

Particle Release:

The filtrate complies with USP <788> test for particulate matter in injections

TOC and Conductivity

Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after a 3 liter flush

Extractables with WFI

Passes NVR test as per USP <661>

Oxidizable Substances

Passes test as per USP <1231>

Ordering Information

Type		Size		Pore Size		Inlet /Outlet		X	X	Sterility		Pack Size	
	Code	Length and EFA	Code		Code		Code				Code		Code
AseptiCap WS (0.45 µm Upstream)	DWSX	1" (250 cm ²)	51	0.2 µm	01	¼" SHB	A			Non-Sterile	1	1	01
		2" (500 cm ²)	52	0.45 µm	02	½" Hose Barb	D			EO Sterile	2		
AseptiCap WS (0.8 µm Upstream)*	DWS5	5" (1000 cm ²)	53			1½" Sanitary Flange	E						
		8" (2000 cm ²)	57			¾" Sanitary Flange	S						
						Quick Connector	J						
						Single Step ½" Hose Barb**	Q						
						Female Luer Lock	U						
						Male Luer Slip***	W						
						3/16" Hose Barb****	N						
						3/8" Hose Barb**	I						
Example		DWSX	53	01		EE		X	X	1		01	

EFA: Effective Filtration Area

Example

*0.8 µm upstream is available with 0.45 µm capsule filters only

**Single step ½" hose barb and 3/8" Hose Barb end connections are not available in 1" capsule filter

***Male luer slip end connection is available as outlet only in 1" capsule filters

****3/16" hose barb end connection is available in:

- 1" and 2" capsule filters as inlet and outlet
- 5" as outlet only